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

AUDIOMETRÍA TONAL LIMINAR

The image shows a screenshot of the Kiverseal audiometry software interface. The interface is divided into two main sections: 'Derecho' (Right) and 'Izquierdo' (Left), each with a frequency-response graph. The central control panel includes buttons for 'Tono HL', ear selection ('R' for Right, 'L' for Left), frequency selection ('1 K Hz'), 'Masking', and volume control. At the bottom, there are buttons for 'Enviar señal', 'Cambio de intensidad', 'Hablar con el paciente', and 'Oír lo que se manda manda al paciente'. The interface also displays 'Paciente' information, 'ACTIVOS' (A, B, R, L, S), and a date/time stamp (04/04/2024 09:17).

Callout boxes and arrows point to the following features:

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
- 1. 7" color graphic TFT display module.
- 2. Built-in thermal printer.
- 3. Integrated microphone module for patient talk over (automatically excluded upon hookup of the external microphone into the appropriate jack connector (26) located on the back panel).
- 4. Channel 1 (dB) and 2 (dB) adjustment knobs: the channel 1 adjustment knob enables channel 1 dB level adjustments and selection of the items on the menu pages. The channel 2 adjustment knob enables channel 2 dB level adjustments and selection of the menu page items. Both adjustment knobs enable sound level controls via the volume control menu page (CONTROLS key 14).
- 5. F1....F6 function keys: enable execution of the function displayed on the bottom

- section of the screen frame, just above the corresponding keys. Said functions will vary according to the type of test being run or the menu that is displayed.
- 6. TEST: when in testing mode, this key enables return to the main selection page.
- 7. ENTER key for accessing the tests from the main menu.
- 8. RIGHT - LEFT inputs for channels 1 and 2: these keys enable selection of the ear being tested for both channels. The selected stimulus signal is thus sent to the transducer corresponding to the selected ear.
- 9. OUTPUT - INPUT for channels 1 and 2: enables individual selection of the type of stimulus and the transducer used for testing.
- 10. ARROW KEYS: used for the various selections provided on the menu pages. For multifrequency audiometric testing they are used to set and enter the frequency setting.
- 11. NORM/REV for channels 1 and 2: these keys enable stimulus transmission to the transducers. Normal mode (output signal present when pressing SIGNAL key 12). Reverse mode (output signal not present, when pressing SIGNAL key 12).
- 12. SIGNAL for channels 1 and 2 (touch keys): in audiometry testing, with the channel in normal mode, these keys enable transmission of the respective channel signals to the transducers for the length of time they are pressed. In Reverse mode, these keys deactivate the transducer output signal for the length of time they are pressed.
- 13. FREQUENCY: these keys enable selection of the signal frequency sent to the transducers.
- 14. CONTROLS: this key provides access to the external signal adjustment and calibration page, for the monitor-headset, for the patient-to-operator mike (e.g. for talk back) and for the operator-to-patient mike (e.g. for talk over).
- 15. T. OVER: this provides direct access to the talk over microphone volume adjustment controls (for operator to patient communication).
- 16. START/STOP: this key enables start/stop of some of the test options (for Autothreshold, SISI, Tone Decay, Bekesy and Speech Audiometry in clinical audiometry testing).
- 17. STORE: this key enables storage and display of the values measured during audiogram testing.
- 18. PRINT: if pressed from the main selection page, this key enables printing of all the tests performed (via thermal printer, or by creating a pdf report), otherwise printout of a single test is possible if the key is pressed directly from the test page (with printout via the built-in thermal printer).

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